

Human GHRL/MTLRP protein (Recombinant) (C-RaFC) (STJP001964)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51738
Gene Symbol	GHRL
UniProt ID	GHRL_HUMAN
Immunogen Region	Gly24-Lys117
Immunogen Sequence	GSSFLSPEHQRVQQRKESKK PPAKLQPRALAGWLRPEDGG QAEGAEDLEVRFNAPFDVG IKLSGVQYQQHSQALGKFLQ DILWEEAKEAPADK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

