

Human Microtubule-associated tau/MAPT protein (Recombinant) (C-hFC) (STJP001970)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4137
Gene Symbol	MAPT
UniProt ID	TAU_HUMAN
Immunogen Region	Gln88-Arg230
Immunogen Sequence	QAAAPHTPEIPEGTTAEEAG IGDTPSLEDEAAGHVTQARM VSKSKDGTGSDDKKAKGADG KTKIATPRGAAPPQQKQGAN ATRIPAKTPAPKTPSSGE PPKSGDRSGYSSPGSPGTPG SRSRTPSLPTPTREPCKVA VVR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

