

## Human TNFSF8/CD30L protein (Recombinant) (N-His) (STJP001971)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">944</a>
Gene Symbol	<a href="#">TNFSF8</a>
UniProt ID	<a href="#">TNFL8_HUMAN</a>
Immunogen Region	Gln63-Asp234
Immunogen Sequence	QRTDSIPNSPDNVPLKGGNC SEDLLCILKRAPFKKSWAYL QVAKHLNKTLSWNKDGILH GVR YQDGNLVIQFPGLYFII CQLQFLVQCPNNSVDLKLEL LINKHIKKQALVTVCESGMQ TKHVYQNLSQLLDYLVQVNT TISVNVDTFQYIDTSTFPLE NVLSIFLYSNSD

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

