

Human RETN protein (Recombinant) (N-hFc-His) (STJP001976)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56729
Gene Symbol	RETN
UniProt ID	RETN_HUMAN
Immunogen Region	Ser17-Pro108
Immunogen Sequence	SSKTLCSMEEAINERIQEVA GSLIFRAISSIGLECQSVTS RGDLATCPRGFAVTGCTCGS ACGSWDVRAETTCHCQCAGM DWTGARCCRVQP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

