

Human IL1RAPL-2 protein (Recombinant) (C-hFc & 6*His) (STJP001982)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26280
Gene Symbol	IL1RAPL2
UniProt ID	IRPL2_HUMAN
Immunogen Region	Thr17-Glu356
Immunogen Sequence	TNLKMVSKRNSVDGCDWSV DLKTYMALAGEPVRVKCALF YSYIRTNYSTAQSTGLRLMW YKNKGDLEPIIFSEVRMSK EEDSIWFHSAEAQDSGFYTC VLRNSTYCMKVSMSLTVAEN ESGLCYNSRIRYLEKSEVTK RKEISCPDMDDFKKSDQEPD VVWYKECKPKMWRSIIQKG NALLIQEVQEEDGGNYTCEL KYEGKLVRRTELKVTALLT DKPPKPLFPMENQPSVIDV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

