

Human FGF-18 protein (Recombinant) (C-His) (STJP001988)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8817
Gene Symbol	FGF18
UniProt ID	FGF18_HUMAN
Immunogen Region	Glu28-Ala207
Immunogen Sequence	EENVDFRIHVENQTRARDDV SRKQLRLYQLYSRTSGKHIQ VLGRRISARGEDGDKYAQLL VETDTFGSQVRIKGGKETEFY LCMNRKGLVKGPDGTSKEC VFIEKVLNNYTALMSAKYS GWYVGFTKKGRPRKGPKTRE NQQDVHFMKRYPKGQPELQK PFKYTTVTKRSRIRPHTPA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------

