

Human Angiopoietin-1/ANG-1/ANGPT1 protein (Recombinant) (C-His) (STJP001989)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	284
Gene Symbol	ANGPT1
UniProt ID	ANGP1_HUMAN
Immunogen Region	Arg277-Phe498
Immunogen Sequence	KREEEKPF R DCADVYQAGFN KSGIYTIYINNMP E PKKVFC NMDVNGGGWTVIQHREDGSL DFQ R GWKEYKMGFGNPSGEY WLGNEFIFAITSQRQYMLRI ELMDWEGNRAYSQYDRFHIG NEKQNYR L Y L KGHTGTAGKQ SSLILHGADFSTKDADNDNC MCKCALMLTGGWWFDACGPS NLNGMFYTAGQNHGKLNLIK WHYFKGPSYSLRSTTMMIRP LDF

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

