

Human Gastrin-releasing peptide/GRP protein (Recombinant) (C-hFc) (STJP002000)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2922
Gene Symbol	GRP
UniProt ID	GRP_HUMAN
Immunogen Region	Arg22-Gly148
Immunogen Sequence	RAVPLPAGGGTVLTMYPRG NHWAVGHLMGKKSTGESSV SERGSLKQQLREYIRWEEAA RNLLGLIEAKENRNHQPPQP KALGNQQPSWDSSEDSSNFKD VGSKGKVGRLSAPGSQREGR NPQLNQQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

