

Mouse IL-6RA/CD126 protein (Recombinant) (C-His) (STJP002005)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16194
Gene Symbol	Il6ra
UniProt ID	IL6RA_MOUSE
Immunogen Region	Leu 20-Glu 357
Immunogen Sequence	LVLGSCRALEVANGTVTSLP GATVTLICPGKEAAGNVTIH WVYSGSQNREWTTTGNTLVL RDVQLSDTGDYLCSLNDHLV GTVPLLVDVPPEEPKLSCFR KNPLVNAICEWRPSSTPSPT TKAVLFAKKINTTNGKSDQK VPCQYSQQLKSFSCQVEILE GDKVYHIVSLCVANSVGSKS SHNEAFHSLKMQPDPPANL VVSAIPGRPRWLKVSQWHPH TWDPSSYLLQFQLRYRPVW

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

