

Human ACE/CD143 protein (Recombinant) (N-his) (STJP002006)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1636
Gene Symbol	ACE
UniProt ID	ACE_HUMAN
Immunogen Region	Leu30-Leu1261
Immunogen Sequence	LDPGLQPGNFSADEAGAQLF AQSYNSSAEQVLFQSVAASW AHDTNITAENARRQEEAALL SQEFAEAWGQKAKELYEPIW QNFTDPQLRRIIGAVRTLGS ANLPLAKRQQYNALLSNMSR IYSTAKVCLPNKTATCWSLD PDLTNILASSRSYAMLLFAW EGWHNAAGIPLKPLYEDFTA LSNEAYKQDGFDTGAYWRS WYNSPTFEDDLEHLYQQLEP LYLNLFHAFVRRALHRRYGD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human ACE/CD143 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 140-180 kDa.