

Human C1QC protein (Recombinant) (N-6His & N-hFc) (STJP002008)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 M filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	714
Gene Symbol	C1QC
UniProt ID	C1QC_HUMAN
Immunogen Region	Asn29-Asp245
Immunogen Sequence	NTGCYIGIPGMPGLPGAPGKD GYDGLPGPKGEPGIPAIPGI RGPKGQKGEPGLPGHPGKNG PMGPPGMPGVPGPMGIPGEP GEEGRYKQKFQSVFTVTRQT HQPPAPNSLIRFNAVLTPQ GDYDTSTGKFTCKVPGLYYF VYHASHTANLCVLLYRSGVK VVTFCGHTSKTNQVNSGGVL LRLQVGEEVWLAVNDYYDMV GIQGSDSVFSGLLPD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human C1QC Protein was determined by SDS-PAGE with Coomassie Blue, showing bands at 45-50 kDa.