

Human TRAIL-R1/CD261 protein (Recombinant) (C-His) (STJP002037)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8797
Gene Symbol	TNFRSF10A
UniProt ID	TR10A_HUMAN
Immunogen Region	Ala109-Asn239
Immunogen Sequence	ATIKLHDQSIGTQQWEHSPL GELCPPGSHRSEHPGACNRC TEGVGYTNASNNLFACLPCT ACKSDEEEERSPCTTTRNTAC QCKPGTFRNDSAEMCRKCS RGCPRGMVKVKDCTPWSIE CVHKESGNGHN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

