

Mouse CD163 protein (Recombinant) (C-His) (STJP002051)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	93671
Gene Symbol	Cd163
UniProt ID	C163A_MOUSE
Immunogen Region	Met1-Thr1045
Immunogen Sequence	MGGHRMVLLGGAGSPGCKRF VHLGFFVAVSSLLSASAVT NAPGEMKKELRLAGGENNCS GRVELKIHDKWGTVCNNGWS MNEVSVVCQQLGCPTSIAL GWANSSAGSGYIWMKVSCT GNESALWDCKHDGWGKHNT HEKDAGVTCSDGSNLEMRLV NSAGHRCLGRVEIKFQGWG TVCDNFSKDHASVICKQLG CGSAISFSGSAKLGAGSGPI WLDDLACNGNESALWDCKH

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

