

## Human Decorin/PG-S2/DCN protein (Recombinant) (C-hFC) (STJP002063)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1634</a>
Gene Symbol	<a href="#">DCN</a>
UniProt ID	<a href="#">PGS2_HUMAN</a>
Immunogen Region	Asp31-Lys359
Immunogen Sequence	DEASGIGPEVPDDRDFEPSL GPVCPFRQCQCHLRVVQCSDL GLDKVPKDLPPDTLLDLQN NKITEIKDGDGFKNLKLNHAL ILVNNKISKVSPGAFTPLVK LERLYLSKNQLKELPEKMPK TLQELRAHENEITKVRKVTF NGLNQMIVIELGTNPLKSSG IENGAFQGMKKLSYIRIADT NITSIPQGLPPSLTELHLDG NKISRVDAAASLKGLNNAKL GLSFNSISAVDNGSLANTP

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

