

## Human Angiopoietin-like 4/ANGPTL4 protein (Recombinant) (N-His) (STJP002072)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">51129</a>
Gene Symbol	<a href="#">ANGPTL4</a>
UniProt ID	<a href="#">ANGL4_HUMAN</a>
Immunogen Region	Pro166-Ser406
Immunogen Sequence	PEMAQPVDPAHNVSRLHRLP RDCQELFQVGERQSGLFIEIQ PQGSPFFLVNCKMTSDGGWT VIQRRHDGSDVFNRPWEAYK AGFGDPHGFEWLGLEKVHSI TGDRNSRLAVQLRDWDGNAE LLQFSVHLGGEDTAYSLQLT APVAGQLGATTVPPSGLSVP FSTWDQDHLRRDKNCAKSL SGGWWFGTCSHSNLNGQYFR SIPQQRQKLKKGIFWKTWRG RYYPLQATTMLIQPMAAEA

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------

