

Human SR-B3/CD36 protein (Recombinant) (N-hFc) (STJP002075)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	948
Gene Symbol	CD36
UniProt ID	CD36_HUMAN
Immunogen Region	Gly30-Asn439
Immunogen Sequence	GDLLIQKTIKKQVLEEGLTI AFKNWVKTGTEVYRQFWIFD VQNPQEVMMNSSNIQVKQRG PYTYRVRFLAKENVTQDAED NTVSFLQPNGAIFEPSLSVG TEADNFTVLNLAVAAASHIY QNQFVQMILNSLINKSKSSM FQVRTLRELLWGYRDPFLSL VPYPVTTTGLFYYPYNNNTAD GVYKVFNGKDNISKVAIIDT YGKRNL SYWESHCDMINGT DAASFPPFVEKSQVLQFFS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

