

Human Fas ligand/APTL/CD95 ligand/CD178 protein (Recombinant) (N-His) (STJP002146)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	356
Gene Symbol	FASLG
UniProt ID	TNFL6_HUMAN
Immunogen Region	Pro134-Leu281
Immunogen Sequence	PSPPEKELRKAHVAHTGKSNRSRSMPLWEDTYGIVLLSGVKYKKGGLVINETGLYFVYS KVYFRGQSCNNLPLSHKVYM RNSKYPQDLVMMEGKMMSYC TTGQMWARSSYLGAVFNLTS ADHLYVNVSELSLVNFEESQ TFFGLYKL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

