

## Human Flt3 ligand/Flt3L protein (Recombinant) (C-hFc) (STJP002148)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2323</a>
Gene Symbol	<a href="#">FLT3LG</a>
UniProt ID	<a href="#">FLT3L_HUMAN</a>
Immunogen Region	Thr27-Pro185
Immunogen Sequence	TQDCSFQHSPISSDFAVKIR ELSDYLLQDYPVTVASNLQD EELCGGLWRLVLAQRWMERL KTVAGSKMQGLLERVNTIEIH FVTKCAFQPPPSCLRFVQTN ISRLLQETSEQLVALKPWIT RQNFSRCLELQCQPDSSSTLP PPWSPRPLEATAPTAPQPP

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

