

Human CRTAM/CD355 (A65V) protein (Recombinant) (C-His) (STJP002331)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Gene ID	56253
Gene Symbol	CRTAM
UniProt ID	CRTAM_HUMAN
Immunogen Region	Ser18-Gly287
Immunogen Sequence	LTNHTETITVEEGQTLTLKC VTSLRKNSSLQWLTPSGFTI FLNEYPLVKNskyQLLHHSANQLSITVPNVTLQDEGVYKC LHYSDSVSTKEVKVIVLATP FKPILEASVIRKQNGEEHVV LMCSTMRSKPPPPQITWLLGN SMEVSGGTLHEFETDGKKCN TTSTLIHTYGKNSTVDCII RHRGLQGRKLVAPFRFEDLV TDEETASDALERNSSLSSQDP QQPTSTSVSTEDSSTSEID

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

