

Mouse NKG2D/CD314 protein (Recombinant) (N-hFc) (STJP002377)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Gene ID	27007
Gene Symbol	Klrk1
UniProt ID	NKG2D_MOUSE
Immunogen Region	Phe90-Val232
Immunogen Sequence	FKETFQPVLCNKEVPVSSRE GYCGPCPNWICHNRNICYQF FNEEKTWNQSQASCLSQNSS LLKIYSKEEQDFLKLKLSYH WMGLVQIPANGSWQWEDGSS LSYNQLTLVEIPKGSCAVYG SSFKAYTEDCANLNTYICMK RAV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

