

Human Tenascin protein (Recombinant) (C-His) (STJP002430)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3371
Gene Symbol	TNC
UniProt ID	TENA_HUMAN
Immunogen Region	Gly23-Ser621
Immunogen Sequence	GVLKKVIRHKRQSGVNATLP EENQPVVFHVNLIKLPVGS QCSVDLESASGEKDLAPPSE PSESFQEHTVDGENQIVFTH RINIPRRACGCAAAPDVKEL LSRLEELENLVSSLREQCTA GAGCCLQPATGRLDTRPFCS GRGNFSTEGCGVCEPGWKG PNCSEPECPGNCHLRGRCID GQCICDDGFTGEDCSQLACP SDCNDQGKCVNGVCICFEGY AGADCSREICPVCSEEHG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

