

Human Cathepsin B protein (Recombinant) (C-His) (STJP002444)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 Cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1508
Gene Symbol	CTSB
UniProt ID	CATB_HUMAN
Immunogen Region	Arg18-Ile339
Immunogen Sequence	RSRPSFHPLSDELVNYVNKR NTTWQAGHNFYNVDMSYLKR LCGTFLGGPKPPQRMFTED LKLPA SFDAREQWPQCPTIK EIRDQGSCGSCWAFGAVEAI SDRICHTNAHVSVEVSAED LLTCCGSMCGDGCNGGYP AE AWNFWTRKGLVSGGLYESHV GCRPYSIPPCEHHVNGSRPP CTGEGDTPKCSKICEPGYSP TYKQDKHYGYNYSVSNSEK DIMAEIYKNGPVEGA FSVY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

