

Human Angiopoietin-like 3/ANGPTL3 protein (Recombinant) (C-His) (STJP002462)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, 0.05% CHAPS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27329
Gene Symbol	ANGPTL3
UniProt ID	ANGL3_HUMAN
Immunogen Region	Ser17-Glu460
Immunogen Sequence	SRIDQDNSSFDSLSPEPKSR FAMLDDVKILANGLLQLGHG LKDFVHKTKGQINDIFQKLN IFDQSFYDLSLQTSEIKEEE KELRRTTYKLQVKNEEVKNM SLELNSKLESLEEKILLQQ KVKYLEEQLTNLIQNQPETP EHPEVTSKTFVEKQDNSIK DLLQTVEDQYQLNQQHSQI KEIENQLRRTSIQEPTEISL SSKPRAPRTPFLQLNEIRN VKHDGIPAECTTIYNRGEH

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

