

## Human CCL18/PARC protein (Recombinant) (No-Tag) (STJP002480)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	Yeast

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6362</a>
Gene Symbol	<a href="#">CCL18</a>
UniProt ID	<a href="#">CCL18_HUMAN</a>
Immunogen Region	Ala21-Ala89
Immunogen Sequence	AQVGTNKELCCLVYTSWQIP QKFIVDYSETSPQCPKPGVI LLTKRGRQICADPNKKWVQK YISDLKLNA

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

