

Human GAS-6 protein (Recombinant) (C-His) (STJP002482)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2621
Gene Symbol	GAS6
UniProt ID	GAS6_HUMAN
Immunogen Region	Ala49-Ala678
Immunogen Sequence	MAPSLSPGPAALRRAPQLLL LLLAAECALAALLPAREATQ FLRPRQRRRAFQVFEEAKQGH LERECVEELCSREEAREVFE NDPETDYFYPRYLDCKINKYG SPYTKNSGFATCVQNLDPQC TPNPCRDRKGTQACQDLGMNF FCLCKAGWGGRLCDKDVNEC SQENGGCLQICHNKPGSFHC SCHSGFELSSDGRTCQDIDE CADSEACGEARCKNLPGSYS CLCDEGFAYSSQEACRDV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

