

Human NCAM-1/CD56 protein (Recombinant) (C-His) (STJP002495)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4684
Gene Symbol	NCAM1
UniProt ID	NCAM1_HUMAN
Immunogen Region	Leu20-Gly718
Immunogen Sequence	LQVDIVPSQGEISVGESKFF LCQVAGDAKDKDISWFSPNG EKLTPNQQRISVVWDDSSS TLTIYNANIDDAgiYKCVVT GEDGSESEATVNVKIFQKLM FKNAPTPQEFREGEDAVIVC DVVSSLPPTIIWKHKGRDVI LKKDVRFIVLSNNYLQIRGI KKTDEGTYRCEGRILARGEI NFKDIQVIVNPPTIARQN IVNATANLGQSVTLVCDAEG FPEPTMSWTKDGEQIEQEE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

