

Aequorea victoria EGFP protein (Recombinant) (C-His) (STJP002540)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Met1-Lys238
Immunogen Sequence	MVSKGEELFTGVVPILVELD GDVNGHKFVS SGEGEDATY GKLT LKFICTTGKLPVPWPT LVTTLTYGVQCFSRYPDHMK QHDFFKSAMPEGYVQERTIF FKDDGNYKTRAEVKFEGDTL VNRIELKGIDFKEDGNILGH KLEYNYNSHNVYIMADKQKN GIKVNFKIRHNIEDGSVQLA DHYQQNTPIGDGPVLLPDNH YLSTQSALS KDPNEKRDHMV LLEFVTAAGITLGMDELYK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

