

Human IgG3 protein (Recombinant) (No-Tag) (STJP002543)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene Symbol	IGHG3 Uniprot ID= IGHG3_HUMAN "
Immunogen Region	Glu99-Lys377
Immunogen Sequence	ELKTPLGDTTHTCPRCPEPK SCDTPPPCPRCPEPKSCDTP PPCPRCPEPKSCDTPPPCPR CPAPELLGGPSVFLFPPKPK DTLMSRTPEVTCVVVDVSH EDPEVQFKWYVDGVEVHNAK TKPREEQYNSTFRVSVLTV LHQDWLNGKEYKCKVSNKAL PAPIEKTISKTKGQPREPQV YTLPPSREEMTKNQVSLTCL VKGFYPSDIAVEWESSGQPE NNYNTTPMLSDGSFFLY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

