

Human IgG4 protein (Recombinant) (C-His) (STJP002544)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene Symbol	IGHG4 Uniprot ID= IGHG4_HUMAN "
Immunogen Region	Glu99-Lys327
Immunogen Sequence	ESKYGPPCPSCPAPEFLGGP SVFLFPPKPKDTLMISRTPE VTCVVVDVSQEDPEVQFNWY VDGVEVHNAKTKPREEQFNSTYRVVSVLTVLHQDWLNGKE YKCKVSNKGLPSSIEKTIK AKGQPREPQVYTLPPSQEEM TKNQVSLTCLVKGFYPSDIA VEWESNGQPENNYKTTTPVL DSDGSFFLYSRLTVDKSRWQ EGNVFSCSVMEALHNHYTQ KLSLSLGLK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

