

Human IgE protein (Recombinant) (C-His) (STJP002548)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene Symbol	IGHE Uniprot ID= IGHE_HUMAN "
Immunogen Region	Ser106-Gly427
Immunogen Sequence	SRDFTPPTVKILQSSCDGGG HFPPTIQLLCLVSGYTPGTI NITWLEDGQVMDVDLSTAST TQEGELASTQSELTLSQKHW LSDRTYTCQVTYQGHTFEDS TKKCADSNPRGVSA YLSRPS PFDLFIRKSP TITCLVVDLA PSKGTVNLTWSRASGKPVNH STRKEEKQRNGTLTVTSTLP VGTRDWIEGETYQCRVTHPH LPRALMRSTTKTSGPRAAPE VYAFATPEWPGSRDKRTLA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

