

Human ErbB-3/HER3 protein (Recombinant) (C-His) (STJP002556)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2065
Gene Symbol	ERBB3
UniProt ID	ERBB3_HUMAN
Immunogen Region	Ser20-Thr643
Immunogen Sequence	SEVGNSQAVCPGTLNGLSVT GDAENQYQTLTKLYERCEVV MGNLEIVLTGHNADLSFLQW IREVTGYVLVAMNEFSTLPL PNLRVVRGTQVYDGKFAIFV MLNYNTNSSHALRQLRLTQL TEILSGGVYIEKNDKLCHMD TIDWRDIVRDRDAEIVVKDN GRSCPPCHEVCKGRCWGPGS EDCQTLTKTICAPQCNGHCF GPNPNQCCHDECAGGCSGPQ DTDCFACRHFNDSGACVPR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human ErbB-3/HER3 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 100 kDa.



Immobilized Recombinant Human ERBB3 at 2 Mu g/mL (100 Mu L/well) can bind Recombinant Human NRG1-beta1 with a linear range of 0. 128-45. 7 ng/mL.