

Human Kallikrein-11/KLK11 protein (Recombinant) (C-His) (STJP002583)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M µm filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11012
Gene Symbol	KLK11
UniProt ID	KLK11_HUMAN
Immunogen Region	Met1-Asn250
Immunogen Sequence	MRILQLILLALATGLVGGET RIIKGFECKPHSQPWQAALF EKTRLLCGATLIAPRWLLTA AHCLKPRYIVHLGQHNLQKE EGCEQTRTATESFPHPGFNN SLPNKDHRNDIMLVKMASPV SITWAVRPLTLSSRCVTAGT SCLISGWGSTSSPQLRLPHT LRCANITIIHQKCENAYPG NITDTMVCASVQEGGKDSCQ GDSGGPLVCNQSLQGIISWG QDPCAITRKPGVYTKVCKY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---