

Human NCR3/NKp30/CD337 protein (Recombinant) (C-hFc & His) (STJP002593)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	259197
Gene Symbol	NCR3
UniProt ID	NCTR3_HUMAN
Immunogen Region	Leu19-Thr138
Immunogen Sequence	LWVSQPPEIRLTLEGSSAFLP CSFNASQGRLAGSVTWFRD EVVPGKEVRNGTPEFRGRLA PLASSRFLHDHQAEHLHVDV RGHDAIYVCRVEVLGLGVG TGNGTRLVVEKEHPQLGAGT

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human NCR3/NKp30/CD337 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 54-60 kDa.



The purity of Human NCR3/NKp300 Protein (Cat. STJP002593) was greater than 95% as determined by SEC-HPLC.



Immobilized recombinant Human B7-H6 at 2 Mu g/mL (100 Mu L/well) can bind recombinant Human NCR3, the EC₅₀ of Human NCR3 is 66.43 ng/mL.