

## Human IFN-lambda 3/IL-28B protein (Recombinant) (C-His) (STJP002601)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES


|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilised from a 0.22 Mu m filtered solution of PBS, 500mM NaCl, pH 7.4.                     |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION


|                    |                                                                                                                                                                                                                     |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">282617</a>                                                                                                                                                                                              |
| Gene Symbol        | <a href="#">IFNL3</a>                                                                                                                                                                                               |
| UniProt ID         | <a href="#">IFNL3_HUMAN</a>                                                                                                                                                                                         |
| Immunogen Region   | Met1-Val196                                                                                                                                                                                                         |
| Immunogen Sequence | MTGDCMPVLVMAAVLTVTG AVPVARLRGALPDARGCHIA QFKSLSPQELQAFKRAKDAL<br>EESLLKDKCKRSRLFPRTW DLRQLQVRERPVALEALAL TLKVLEATADTDPALGDVLD<br>QPLHTLHHLSQLRACIQPQ PTAGPRTRGRLLHHLHRLQE APKKESPGCLEASVTFNLFR<br>LLTRDLNCVASGDL CV |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human IFN-lambda 3/IL-28B Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 23 kDa.



Immobilized Recombinant Human IL28B at 2 Mu g/mL (100 Mu L/well) can bind Recombinant Human IL10RB with a linear range of 0. 25-1 Mu g/mL.