

Human Growth hormone receptor/GHR protein (Recombinant) (C-hFc & His) (STJP002605)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES



Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2690
Gene Symbol	GHR
UniProt ID	GHR_HUMAN
Immunogen Region	Ala27-Tyr264
Immunogen Sequence	AILSRAPWSLQSVNPLKTN SSKEPKFTKCRSPERETFSC HWTDEVHHGTKNLGPIQLFY TRRNTQEWTEWKECPDYVS AGENSCYFNSSFTSIWIPYC IKLTSSGGTVDEKCFVDEI VQPDPIALNWTLLNVSLTG IHADIQVRWEAPRNADIQKG WMVLEYELQYKEVNETKWKM MDPILTTSVPVYSLKVDKEY EVRVRSKQRNSGNYGEFSEV LYVTLPQMSQFTCEEDFY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human Growth hormone receptor/GHR Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 70-85 kDa.

Immobilized Human GH at 2 Mu g/mL (100 Mu L/well) can bind Human GHR with a linear range of 0.1-8.58 ng/mL.