

Human IL1RL1/ST2 protein (Recombinant) (C-His) (STJP002609)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9173
Gene Symbol	IL1RL1
UniProt ID	ILRL1_HUMAN
Immunogen Region	Lys19-Phe328
Immunogen Sequence	KFSKQSWGLENEALIVRCPR QGKPSYTVDWYYSQTNKSIP TQERNRVFASGQLLKFLPAA VADSGIYTCIVRSPTFNRTG YANVTIYKKQSDCNVPDYL M YSTVSGSEKNSKIYCPIDL YNWTAPLEWFKNCQALQGS R YRAHKSFLVIDNVMTE DAGD YTCKFIHNENGANYSVTATR SFTVKDEQGFSLFPVIGAPA QNEIKEVEIGKANLTC SAC FGKGTQFLAAVLWQLNGTK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human IL1RL1/ST2 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 45-50kDa.



Immobilized Recombinant Human IL33, No Tag at 10 µg/mL (100 µL/well) can bind Recombinant Human ST2, His Tag with a linear range of 1.41-5.64 µg/mL.