

Human TNFRSF9/4-1BB/CD137 protein (Recombinant) (C-hFc & His) (STJP002612)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3604
Gene Symbol	TNFRSF9
UniProt ID	TNR9_HUMAN
Immunogen Region	Leu 24-Gln 186
Immunogen Sequence	LQDPCSNCPAGTFCDNRRNQ ICSPCPPNSFSSAGGQRTCD ICRCQCKGVFTRKECSSTSN AECDCPTPGFHC LGAGCSMCE QDCKQGQELTKKGCKDCCFG TFNDQKRGICRPWTNCSLDG KSVLVNGTKERDVVCGPSPA DLSPGASSVTPPAPAREPGH SPQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human TNFRSF9/4-1BB/CD137 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 60-66 kDa.



Immobilized Human TNFRSF9 at 2 Mu g/mL (100 Mu L/well) can bind Human TNFRSF9 with a linear range of 0.1-12.9 ng/mL.