

Human SLAMF1/CD15 protein (Recombinant) (C-hFC) (STJP002630)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6504
Gene Symbol	SLAMF1
UniProt ID	SLAF1_HUMAN
Immunogen Region	Ala21-Pro237
Immunogen Sequence	ASYGTGGRMMNCPKILRLG SKVLLPLTYERINKSMNKSI HIVVTMAKSLNSVENKIVS LDPSEAGPPRYLGDYKFYL ENLTGIRESRKEDEGWYLM TLEKNVSVQRFCLQLRLYEQ VSTPEIKVLNKTQENGCTL ILGCTVEKGDHVAYSWSEKA GTHPLNPANSSHLLSLTLGP QHADNIYICTVSNPISNNSQ TFSWPWPGCRTDPSETKP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---