

Human IFNA6/IFN-alpha 6/LeIF K protein (Recombinant) (C-6xHis) (STJP002634)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3443
Gene Symbol	IFNA6
UniProt ID	IFNA6_HUMAN
Immunogen Region	Ser21-Glu189
Immunogen Sequence	SLDCDLPQTHSLGHRRTMML LAQMRRISLFSCLKDRHDFR FPQEEFDGNQFQKAEAISVL HEVIQQTFNLFSTKDSSVAW DERLLDKLYTELYQLNDLE ACVMQEVWVGGTPLMNEDSI LAVRKYFQRITLYLTEKKYS PCAWEVVRAEIMRSFSSSRN LQERLRKE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---