

Human BST-1/CD157 protein (Recombinant) (C-His) (STJP002635)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	683
Gene Symbol	BST1
UniProt ID	BST1_HUMAN
Immunogen Region	Gly29-Lys292
Immunogen Sequence	GARARWRGEGTSAHLRDIFL GRCAEYRALLSPEQRNKNCT AIWEAFKVALDKDPCSVLPS DYDLFINLSRHSIPRDKSLF WENSHLLVNSFADNTRRFMP LSDVLYGRVADFLSWCRQKN DSGLDYQSCPTSEDCENNPV DSFWKRASIYQSKDSSGVH VMLNGSEPTGAYPIKGFFAD YEIPNLQKEKITRIEIVVMH EIGGPNVESC GEGSMKVLEK RLKDMGFQYSCINDYRPVK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---