

Mouse CD47/IAP/OA3 protein (Recombinant) (C-hFc) (STJP002655)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	16423
Gene Symbol	Cd47
UniProt ID	CD47_MOUSE
Immunogen Region	Gln19-Lys140
Immunogen Sequence	QLLFSNVNSIEFTSCNETVV IPCIVRNVEAQSTEEMFVKW KLNKSYIFIYDGNKNSTTTD QNFTSAKISVSDLINGIASL KMDKRDAMVGNYTCEVTELS REGKTVIELKNRTVSWFSPN EK
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.