

Human IL-17F protein (Recombinant) (N-His) (STJP002672)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	112744
Gene Symbol	IL17F
UniProt ID	IL17F_HUMAN
Immunogen Region	Arg31-Gln163
Immunogen Sequence	RKIPKVGHTFFQKPESCPPV PGGSMKLDIGIINENQRVSM SRNIESRSTSPWNYVTWDP NRYPSEVVQAQCRNLGCINA QGKEDISMNSVPIQQETLVV RRRKHQGCVSFQLEKVLVTV GCTCVTPVIHRVQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---