

## Human Jagged1/JAG1/CD339 protein (Recombinant) (C-hFc) (STJP002673)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">182</a>
Gene Symbol	<a href="#">JAG1</a>
UniProt ID	<a href="#">JAG1_HUMAN</a>
Immunogen Region	Trp896-Arg1065
Immunogen Sequence	WCGPRPCLLHKHSECPSPGQ SCIPILDDQCFVHPCTGVGE CRSSSLQPVKTKCTSDSYQQ DNCANITFTFNKEMMSPGLT TEHICSELRNLNLIKNSAE YSIYIACEPSPSANNEIHVA ISAEDIRDDGNPIKEITDKI IDLVSKRDGNSLIAAVAEV RVQRRPLKNR

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---