

Human CEACAM6/CD66c protein (Recombinant) (C-His) (STJP002676)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4680
Gene Symbol	CEACAM6
UniProt ID	CEAM6_HUMAN
Immunogen Region	Lys35-Gly320
Immunogen Sequence	KLTIESTPFNVAEGKEVLLL AHNLPQNRIGYSWYKGERVD GNSLIVGYVIGTQQATPGPA YSGRETIYPNASLLIQNVQT NDTGFYTLQVIKSDLVNEEA TGQFHVYPELPKPSISSNNS NPVEDKDAVAFTCEPEVQNT TYLWWVNGQSLPVSPRLQLS NGNMTLTLLSVKRNDAGSYE CEIQNPASANRSDPVTNLNL YGPDGPTISPSKANYRPGEN LNLSCAASNPPAQYSWFI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human CEACAM6/CD66c Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 45-65 kDa.



Immobilized Recombinant Human CEACAM6 at 2 Mu g/mL (100 Mu L/well) can bind Recombinant Human CEACAM8 with a linear range of 0. 1-0. 4 Mu g/mL.