

Human Ephrin-A3/EFNA3 protein (Recombinant) (C-hFc & His) (STJP002684)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1944
Gene Symbol	EFNA3
UniProt ID	EFNA3_HUMAN
Immunogen Region	Met 1-Ser 213
Immunogen Sequence	MAAAPLLLLLLVPVPLLPL LAQGGGALGNRHAVYWNSS NQHLRREGYTVQVNVNDYLD IYCPHYNSSGVGPGAGPGPG GGAEQYVLYMVSRRNGYRTCN ASQGFKRWECNRPHAPHSP KFSEKFQRYSAFSLGYEFHA GHEYIYISTPTHNLHWKCLR MKVFCVCCASTSHSGEKPVT LPQFTMGPNVKINVLEDFEG ENPQVPKLEKSI

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human Ephrin-A3/EFNA3 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 68 kDa.



Immobilized Mouse EPHA6 at 1 Mu g/mL (100 Mu L/well) can bind Human EphrinA3 with a linear range of 0. 3-113. 9ng/mL.



Immobilized Mouse EPHA6 at 1 Mu g/mL (100 Mu L/well) can bind Human EphrinA3 with a linear range of 0. 3-59 ng/mL.