

Human TEM8/ANTXR1 protein (Recombinant) (C-His) (STJP002702)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84168
Gene Symbol	ANTXR1
UniProt ID	ANTR1_HUMAN
Immunogen Region	Glu33-Lys321
Immunogen Sequence	EDGGPACYGGFDLYFILDKS GSVLHHWNEIYYFVEQLAHK FISPQLRMSFIVFSTRGTTL MKLTEDREQIRQGLEELQKV LPGGDTYMHEGFERASEQIY YENRQGGYRTASVIIALTDGE LHEDLFFYSEREANRSRDLG AIVYCVGVKDFNETQLARIA DSK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---