

Human Flt4 ligand/VEGF-C protein (Recombinant) (C-His) (STJP002723)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7424
Gene Symbol	VEGFC
UniProt ID	VEGFC_HUMAN
Immunogen Region	Thr103-Arg227
Immunogen Sequence	TEETIKFAAAHYNTEILKSI DNEWRTQCMPREVCIDVGK EFGVATNTFFKPPCVSVYRC GGCCNSEGLQCMNTSTSYLS KTLFEITVPLSQGPKPVTIS FANHTSCRCMSKLDVYRQVH SIIR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human Flt4 ligand/VEGF-C Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 19-25 kDa..



Immobilized Recombinant Human VEGF-C, His Tag at 0.5 Mu g/mL (100 Mu L/well) can bind Recombinant Human VEGFR3, Fc Tag with a linear range of 3.92-15.70 ng/mL.