

Human FGL1/HP-041 protein (Recombinant) (C-His) (STJP002744)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Supplied as a 0.22 Mu m filtered solution in 20mM Tris, 300mM NaCl, 10% Glycerol, pH 8.0
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2267
Gene Symbol	FGL1
UniProt ID	FGL1_HUMAN
Immunogen Region	Met1-Ile312
Immunogen Sequence	MAKVFSFILVTTALTMGREI SALEDCAQEQMRLRAQVRLL ETRVKQQQVKIKQLQENEV QFLDKGDENTVIDLGSKRQY ADCSEIFNDGYKLSGFYKIK PLQSPAEFVSVYCDMSDGGGW TVIQRSDGSENFNRGWKDY ENGFGNFVQKHGEYWLGNKN LHFLTQEDYTLKIDLADFE KNSRYAQYKNFKVGDEKNFY ELNIGEYSGTAGDSLGNFH PEVQWWASHQRMKFSTWDR

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human FGL1/HP-041 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 36 kDa.



Immobilized recombinant Human FGL1 at 5 Mu g/mL (100 Mu L/well) can bind recombinant Human LAG3 with a linear range of 0. 5-21. 6ng/mL.



Immobilized Recombinant Human FGL1 at 1 Mu g/mL (100 Mu L/well) can bind FGL1 Rabbit mAb with a linear range of 1-27. 5 ng/mL.