

## Human DPP4/CD26 protein (Recombinant) (C-His) (STJP002746)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Formulation         | Supplied as a 0.22 µg/ml filtered solution in PBS, pH 7.4.                                     |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                   |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">1803</a>                                                                                                                                                                                                                                              |
| Gene Symbol        | <a href="#">DPP4</a>                                                                                                                                                                                                                                              |
| UniProt ID         | <a href="#">DPP4_HUMAN</a>                                                                                                                                                                                                                                        |
| Immunogen Region   | Ser39-Pro766                                                                                                                                                                                                                                                      |
| Immunogen Sequence | SRKTYTLTDYLNKNTYRLKLY SLRWISDHEYLYKQENNILV FNAEYGNSSVFLENSTFDEF<br>GHSINDYSISPDGQFILLEY NYVKQWRHSYASYDIYDLN KRQLITEERIPNNTQWVTWS<br>PVGHKLAYVWNNDIYVKIEP NLPSYRITWTGKEDIYNGI TDWVYEEEVFSAYSALWWSP<br>NGTFLAYAQFNDTEVPLIEY SFYSDESLQYPKTVRVPYPK AGAVNPTVKFFVNTDSLS |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|