

Human Fc-gamma RIIa/CD32a (H167R) protein (Recombinant) (C-hFc & His) (STJP002756)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2212
Gene Symbol	FCGR2A
UniProt ID	FCG2A_HUMAN
Immunogen Region	Ala36-Ile218
Immunogen Sequence	AAPPKAVLKLEPPWINVLQE DSVTLTCQGARSPESDSIQW FHNGNLIPHTHTQPSYRFKAN NNDSGEYTCQTGQTSLSDPV HLTVLSEWLVLTQPHLEFQE GETIMLRCHSWKDKPLVKVT FFQNGKSQKFSHLDPTFSIP QANHSHSGDYHCTGNIGYTL FSSKPVRI TVQVPSMGSSSP MGI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human Fc-gamma RIIa/CD32a (H167R) Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 60kDa.



Immobilized Human FCGR2A/CD32a (R167) at 1 Mu g/mL (100 Mu L/well) can bind FCGR2A Rabbit pAb with a linear range of 0.977-175.35 ng/mL.